

aquapick[®]

www.aquapick.com

| **Aquapick Co., Ltd.** |

20F, SK V1 Center, 171, Gasan digital 1-ro, Geumcheon-gu, Seoul, Korea | E-mail. info@aquapick.com

Customer Service: +82-1588-2809

20th
Anniversary
SINCE 2001

aquapick[®]

Officially recommended products by
Korean Dental Association

2023 AQUAPICK CATALOGUE




*The Leader
in oral care solutions*





aquapick®


Company Profile




Lobby



Book Cafe



Office Area



Office Area



Product Demonstration Room



Product Display Room



VIP Reception Room




Conference Hall



Office Area



Island Kitchen



Cafeteria



Aquapick Award and Certification Display



Message from CEO

Although improved living standards have enhanced awareness about the care of the teeth and gums, proper care is still not easy to implement in everyday life. This is primarily due to the presence of blind spots where brushing alone does not offer a good solution. Even brushing three times a day still fails to completely prevent small and large problems with the gums and teeth effectively. These blind spots are areas where a brush cannot reach. Pieces of food left around these areas become hotbeds for growing plaque and accelerating the formation of tartar, presenting a threat to the health of the teeth and gums.

Aquapick, which specializes in developing and manufacturing oral care systems, has made efforts to protect and improve your oral health by offering simple and easy solutions to oral care, which requires more than brushing. Thanks to these efforts, Aquapick has won the Presidential Citation which was given by the President of Korea. This is given to the best of best companies that have contributed to social and economic development. Aquapick has taken the lead in the oral care market and promises to continue to afford customers greater satisfaction with better quality and service through intensive investment in R&D and continuous product innovation.

Aquapick is not satisfied with today's achievements in establishing itself as a leading company in the global oral care market and will continue to make every effort to fulfill its social responsibilities by expanding the base of the oral care system and developing various products, enabling more people to take care of their oral health with ease.

Furthermore, Aquapick will pursue diverse activities for the enhancement of customer satisfaction including fulfillment of thorough warranty services. Aquapick has grown with your continued encouragement. We are thankful for your continued faith and trust and we will press forward in our efforts to meet and exceed your expectations with the highest quality products and services.

Thank you.

CEO of Aquapick Co., Ltd.

Kevin Lee





aquapick®

The products that we are selling are not at their highest quality. It is because we are making better products at the moment.

We are not the 'yes man.'

We are the 'no man.'

We always say 'no.'

We will always say "no" until you say "yes."


Fulfilling customer satisfaction, we will pay back our gratitude with greater customer service.






GOOD HABIT

Correct habits for strong teeth and healthy gums



Removes food debris while
massaging gums gently



Prevents peri-implantitis



Ideal for orthodontic patients



Takes care of pregnant
women's teeth (before & after)



Prevention of
kids' tooth decay



Remove bad breath


Must-have for
your oral health
Aquapick!

Aquapick makes teeth and gums healthy and strong,
Aquapick makes your teeth and gums healthy and strong with best oral care products.



HEALTH

Oral health is the beginning of general health.



→ If you clean your teeth carelessly,



Poor oral hygiene accelerates aging process and reduces cognitive ability.

- Journal of American gerontology (Denmark Research Team) -



Greater risks for pneumonia, asthma, and laryngopharyngitis.

- February, 2002, Issue of the British Journal on the Periodontal Membrane-



Greater risk (up to 2.8 times greater than healthy people) for cardiac arrest.

- International Dental Research Association -



Greater risk (up to 7.8 times greater than healthy people) for premature birth.

- International Dental Research Association -



Infection or inflammation in the roots of teeth, may increase the risk of diabetes by lowering insulin levels and cause difficulty in controlling blood sugar levels.

- Taiwanese Dental Health Association-

RESEARCH


Aquapick is an oral health care device with clinically verified results.

Comparison of interdental plaque using Aquapick

Clinical Institution: The Korean Academy of Oral Health
(done by Paired-test)
Researcher: Dr. Bae
(School of dentistry, Seoul National University)

Aquapick is an oral health care device with clinically verified results.

Interdental plaque index was drastically reduced by using Aquapick.




Comparison between Aquapick, dental floss and interdental brush.

Clinical Institution: The Korean Academy of Oral Health
(done by one-way analysis of variance)
Standard deviation 0.67~1.21
Researcher: Dr. Bae
(School of dentistry, Seoul National University)

Aquapick is an oral health care device with high functionality, convenience and satisfaction.

Aquapick has outstanding functionality and satisfaction after use. It is also very convenient when you clean your teeth that are difficult for dental floss or interdental brush to reach



Aquapick is helpful to prevent gum disease.

aquapick





aquapick®

Aquapick is a daily **must-have**
device for modern people.



| AQ-230



| AQ-120



| AQ-350

NEW

Aquapick Oral Irrigator AQ-350



Officially recommended
by the Korean Dental Association

Product Name: Aquapick Oral Irrigator
Model: AQ-350
Color: Pink / Gray
Components: Main Body, Water Tank, Jet Tip 4EA,
Jet Tip Holder 4EA
Voltage: 100-230V~, 50/60Hz
Power Consumption: 24W
Dimension: 150x90x210mm (WxDxH)
Weight: 710g
Water Tank Capacity: 600ml
Pulsations: 2,400 pulses/minute

Jet Tip Color: ● ● ● ●



**| Pulsating stream of water
at 2,400 pulses/minute**

Effective on removal of food
between tooth & gums



**| Convenient handle
mounting with
ergonomic design**

Effective on removal of food
between tooth & gums



| AUTO START

If you pull up the handle, it
**automatically starts
functioning** by activating
the sensor.
(The world's first development)



| 600ml water tank

Sufficiently usable volume



NEW

Aquapick Oral Irrigator AQ-320




Product Name: Aquapick Oral Irrigator
Model: AQ-320
Components: Main Body, Water Tank,
Jet Tip 4EA, UV Sterilizer
Voltage: 100-240V~, 50/60Hz
Power Consumption: 18W
Dimension: 175x95x190mm (WxDxH)
Weight: 675g
Water Tack Capacity: 800ml
Pulsations: 1,650 pulses/minute

Jet Tip Color: ● ● ● ●

**Delicate water |
pressure control**

12-step water pressure
control dial enables
delicate adjustment



**Ergonomically |
designed handle**

The magnet in the handle
makes it convenient to place
the handle on the body



800ml water tank |

Sufficiently usable volume



| UV Sterilizer

Sterilizes residual bacteria
and virus without using
any chemicals through
ultraviolet(UV) sterilization



| Jet tip cabinet

Convenient and hygienic
storage of 4 jet tips

Aquapick Oral Irrigator AQ-310



Product Name: Aquapick Oral Irrigator
Model: AQ-310
Components: Main Body, Water Tank,
Jet Tip 4EA, Special Tips 3EA
(Tongue scrapper,
Orthodontic tip, Pocket tip)
Voltage: 100-240V~, 50/60Hz
Power Consumption : 12V
Dimension: 140x115x200mm (WxDxH)
Weight: 630g
Water Tack Capacity: 600ml
Pulsations: 1,600 pulses/minute
2 Cleaning Modes: Clean / Massage

Jet Tip Color: ● ● ● ●

Sophisticated water pressure controller

Men and women of all ages
can use it, thanks to
10 steps of sophisticated
water pressure controller.



Jet tip case |

Easily store jet tips
with jet tip case.



600ml water tank |

Sufficiently usable volume



| 3 special tips

Three special tips for
perfect oral care allow
three-dimensional mouthwash



| 4 different colors of clear jet tip

All family members can use it



Aquapick Oral Irrigator AQ-300



Product Name: Aquapick Oral Irrigator

Model: AQ-300

Components: Main Body, Water Tank, Jet Tip 4EA,
Jet Tip Holder 4EA

Voltage: 100-230V~, 50/60Hz

Power Consumption: 24W

Dimension: 150x90x210mm (WxDxH)

Weight: 710g

Water Tank Capacity: 600ml

Pulsations: 2,200 pulses/minute

Jet Tip Color: ● ● ● ●



| Pulsating stream of water at 2,200 pulses/minute

Effective on removal of food between
teeth & gums and massage on the gums



| Water drop shaped tip holders

Sanitarily store jet tips
(can be used as tooth brush holders)



| 4 different colors of clear jet tip

All family members can use it



| Easy control of water pressure

Men and women of all ages can use it
due to sophisticated water pressure control



NEW

Aquapick Cordless Oral Irrigator AQ-250



Product Name: Aquapick Cordless Oral Irrigator

Model: AQ-250

Components: Main Body, Water Tank, Charger, Jet Tip 2EA

Voltage: 5V

Power Consumption: 2W

Dimension: 68x89x288mm (WxDxH)

Weight: 315g

Water Tank capacity: 200ml

Pulsations: 1,600 pulses/min

3 Cleaning Modes: Clean/Soft/Massage

Jet Tip Color: ● ●




• **Strong pulsating stream of water at 1,600 pulses/minute**

Effective on washing away bacteria and food residues in the hidden parts of the mouth together with massage on the gum

• **User-friendly function**

If there is no water in the water tank during use, the power is automatically turned off, so there is no need to press the power button separately to turn it off



• **3 Mode Options**

- Clean** : General cleaning mode
- Soft** : Soft cleaning mode for sensitive teeth
- Massage** : Gum massaging mode

Provide a suitable water pressure for every user, no matter sensitive teeth or not



Aquapick Portable Oral Irrigator AQ-240




Product Name: Aquapick Portable
Oral Irrigator
Model: AQ-240
Components: Main Body, Water Tank, Pouch
Jet Tip 4EA, Charger Cable
Voltage: 5V
Power Consumption: 3W
Dimension: 58x103x132mm (WxDxH)
Weight: 422g
Water Tank capacity: 450ml
Pulsations: 1,400 pulses/minute
2 Cleaning Modes: Clean / Massage

Jet Tip Color: ● ● ● ●



| Pouch for Traveling

Pouch for storing and
protecting the product
when traveling



| Large Capacity Water Tank

Sufficiently use the product
with 450ml water
tank capacity



| Water Pressure Control

Adjust the water pressure in 5 levels
by moving the water pressure control
button up and down



| Portable Size

Neat and small size,
suitable for carrying around
when traveling



Aquapick Cordless Oral Irrigator


AQ-230



Officially recommended
by the Korean Dental Association

Product Name: Aquapick Cordless Oral Irrigator
Model: AQ-230
Components: Main Body, Water Tank, Charger,
Jet Tip 2EA
Voltage: 100-240V~, 50/60Hz
Power Consumption: 2W
Dimension: 60x85x275mm (WxDxH)
Weight: 280g
Water Tank capacity: 200ml
Pulsations: 1,400 pulses/minute
3 Cleaning Modes: Clean / Soft / Massage

Jet Tip Color: ● ●



- Ergonomic jet tips for optimal cleaning

Optimal oral cleaning effect

0.6mm powerful ultra-fine stream

1,400 pulsating streams per minute



- 3 cleaning modes

Clean/Soft/Massage

Available for
customized management



- Conveniently Portable

Optimal size for
carrying around






Aquapick Cordless Oral Irrigator AQ-220



Product Name: Aquapick Cordless Oral Irrigator
Model: AQ-220
Components: Main Body, Water Tank, Charger,
Jet Tip 2EA
Voltage: 110-240V~, 50/60Hz
Power Consumption: 0.7-1.7W
Dimension: 60x75x275mm (WxDxH)
Weight: 260g
Water Tank Capacity: 150ml
Pulsations: 1,400 pulses/minute
3 Cleaning Modes: Normal / Soft / Pulse


Jet Tip Color: ● ●



- 3 Cleaning mode enables customized management
 - Normal** : General cleaning mode
 - Soft** : Soft cleaning mode for sensitive teeth
 - Pulse** : Gum massaging mode
 - CH** : Charging Indicator



- Upgraded charging method
 - Non-contact charging method!**
It will start charging just by placing on the spot.
Easy & safe!



Takes 24 hours
to fully charge



The battery can last up to 2 weeks
when fully charged



NEW

Aquapick Cordless Oral Irrigator AQ-208



Product Name: Aquapick Cordless Oral Irrigator
Model: AQ-208
Color: Pink/White
Components: Main Body, Water Tank, USB Cable,
Pouch, Jet Tip 2EA
Voltage: 5V
Power Consumption: 8W
Dimension: Folded - 30x60x128mm (WxDxH)
Spreaded - 30x60x200mm (WxDxH)
Weight: 200g
Water Tank Capacity: 140ml
Pulsations: 1,600 pulses/minute
3 Cleaning Modes: Clean/Soft/Massage

Jet Tip Color: ● ●



- **Hygienic care**

The wide inlet design and the complete separation of the water tank make it easy to clean




- **Jet tip tray storage bin on the body**

You can conveniently and hygienically store the jet tip



- **Compact size**

Water tank specifically designed to save space
Easy to carry anytime, anywhere





NEW

Aquapick Portable Oral Irrigator AQ-206



Product Name: Aquapick Portable Oral Irrigator
Model: AQ-206
Components: Main Body, Water Tank, USB Cable,
Pouch, Jet Tip 2EA
Voltage: 5V
Power Consumption: 4.5W
Dimension: Folded - 45x65x135mm (WxDxH)
Spreaded - 45x65x205mm (WxDxH)
Weight: 240g
Water Tank Capacity: 180ml
Pulsations: 1,400 pulses/minute
4 Cleaning Modes: Clean/Soft/Massage/Custom

Jet Tip Color: ● ●



• 4 different types of water stream

- Clean** : Clean teeth and gums thoroughly
- Soft** : Clean teeth and gums softly
- Massage** : Massag the gums freshly
- Custom CH** : Choose the water pressure freely



• Jet tip tray storage bin on the body

You can conveniently and hygienically store the jet tip.



• Hygienic care

The wide inlet design and the complete separation of the water tank make it easy to clean.



Aquapick Cordless Oral Irrigator AQ-205



Product Name: Aquapick Cordless Oral Irrigator
Model: AQ-205
Color: White / Pink / Mint
Components: Main Body, Water Tank,
Charging Cable, Jet Tip 2EA
Voltage: 3.7V
Power Consumption: 5W
Dimension: 40.8x58x281.5mm (WxDxH)
Weight: 369g
Water Tank capacity: 170ml
Pulsations: 1,700 pulses/minute
3 Cleaning Modes: Clean / Soft / Massage

Jet Tip Color: ● ●



• Easy to Carry

There is a **jet tip hook** inside the water tank at the bottom of the main body, so it can be easily carried.



• Convenient magnetic charging

USB type cable and **magnetic charging method** make it more **convenient** and **safer** to use than general charging cables.



• Large capacity battery

It can be used for up to **30 days** without charging when it is fully charged for **4 hours**.





Aquapick Sonic Electric Toothbrush AQ-120




Product Name: Aquapick Sonic Electric Toothbrush
Model: AQ-120
Components: Handle, Charger,
Brush Head 2EA
Voltage: 100~240V / 50~60Hz
Power consumption: 0.7~1.5W
Dimension: 30x30x250mm (WxDxH)
Weight: 120g
Vibration: 31,000 pulses/minute
3 Brushing Modes: Clean / Soft / Massage

Brush Head Color: ● ●



- **Portable size**

Convenient to carry around because of the **small size and light weight**



- **Smart Interval Timer**

Smart tooth brush which pauses in **every 30 seconds** It recommends proper brushing habit for children



- **3 different brushing mode**

Clean : General mode for brushing teeth
Soft : Washes sensitive spots or teeth
Massage: Massages gums softly



- **Non-contact charging method**

Designed to be **safe with waterproof** level of IPX 7



NEW

Aquapick Sonic Electric Toothbrush AQ-102




Product Name: Aquapick Sonic
Electric Toothbrush
Model: AQ-102
Color: Midnight Black / Noble white /
Pink Blossom / Mystic Red
Componets: Handle, Charger, Brush Head 2EA
Voltage: 3.7V
Power Consumption : 2W
Dimension: 29.5x29.5x245mm (WxDxH)
Weight: 105g
Vibration: 36,000 pulses/minute
2 Brushing Modes: Clean / Sensitive

Brush Head Color: ● ●



- Sonic vibration bubble with more than 30,000 times/ minute

Gets rid of food debris which is hidden here and there by using sonic waves.



- Smart Interval Timer

Smart tooth brush which pauses in every 30 seconds It recommends proper brushing habit for children



- USB charging type

Convenient to carry around due to small size and light weight



NEW

Aquapick Sonic Electric Toothbrush AQ-101




Product Name: Aquapick Sonic
Electric Toothbrush
Model: AQ-101
Componets: Handle, Charger, Brush Head 2EA
Voltage: 3.7V
Power Consumption : 1.85W
Dimension: 21x21x224mm (WxDxH)
Weight: 88g
Vibration: 40,000 pulses/minute
3 Brushing Modes: Clean/Soft/Massage

Brush Head Color: ● ●



- Dupont company's (U.S.) brush

Minimize damage on teeth enamel by using Dupont company's brush which has the highest quality in the world.



- Smart Interval Timer

Smart tooth brush which pauses in every 30 seconds It recommends proper brushing habit for children



- USB charging type

Small sized and light weighted



Aquapick Toothbrush Charcoal & Silver



• Charcoal Toothbrush


Aquapick Toothbrush Charcoal is an advanced toothbrush containing a charcoal component and having a function of emitting far-infrared rays and removing the odor with a strong adsorption.



4EA



SUPER SOFT
Containing charcoal ingredient in the bristle



SUPER SOFT Fine bristles effectively eliminates microscopic debris and plaque between teeth and gums because one end of the elastic and smooth of the bristle of 0.18mm is processed to a fine cylindrical shape of 0.03mm or less, special functionality.



Silver

• Silver Toothbrush

Aquapick Toothbrush Silver is an advanced toothbrush containing silver component of 99.9% purity and having its strong antibacterial effect proven.



12EA



SUPER SOFT
Containing silver ingredient in the bristle



Aquapick Toothbrush Sterilizer AQ-201



Product Name: Aquapick Toothbrush Sterilizer
Model: AQ-201
Electrical Currency: AC220V / 60Hz
Voltage: 5W
Dimension: 240x62x76mm
Cord Length: 1.3m
Weight: 270g



Heating part is protected by **special Kapton tape**, which is effective for preventing electric shock and blocking electric wave.

- **Semi-permanently usable without changing the lamp**
Easily sterilize toothbrush which is used everyday.
- **24 hours drying & sterilizing germs**
The PTC heater, which uses infrared radiation, sterilize for 24 hours. This keeps toothbrush clean at all times.

Aquapick Portable Toothbrush Sterilizer AQ-203



Product Name: Aquapick Portable Toothbrush Sterilizer
Model: AQ-203
Longevity of Lamp: approximately 20,000 hours
Dimension: 210x50x25mm
Weight: 100g (excluding the battery)



CCGL UV Lamp



Usable for 4~6 month



Wall holder



Convenient cap



ventilation

Aquapick Dental Mask

AQ-50



Product Name: Aquapick Dental Mask
Model: AQ-50
Type: For operation use, For general use
Amount: 50EA (1BOX)
Function: Diagnosis, Infection prevention
Storage: Store at room temperature
(1~30℃)
Dimension: 175x95mm
Registration Number: 5008



BEF

Effective for blocking **over 99.9%** of bacteria

PEF


Effective for blocking **over 99.9%** of particle

FDA

100% non-woven fabric which is harmless



- **FDA** certified non-asbestos material with excellent durability
- | 100% Polypropylene of safe material and passed fire retardant test



- Hypoallergenic, prevents trouble on skin
- | Reduces risk of causing skin trouble because it **doesn't contain any latex, fiberglass, dye, or fluorescent material**



- Fresh air-permeability, soft elasticity
- | **Double wire** provides comfortable wearing sensation
- | Outstanding air permeability enables comfortable breathing



- Triple filter, excellent for blood permeate resistance
- | **Blocks bacteria & dust, prevent respiratory disease.**
- | Filtering efficiency was tested through laboratories in the US
- | * Efficiency certified by the US Nelson Laboratories

Aquapick Propolis Toothpaste

AQ-32



Ministry of Food and
Drug Safety

Officially registered product



Product Name: Aquapick Propolis Toothpaste
Model: AQ-32
Weight: 130g
Storage method: Room temperature storage
(1~30℃)
How to use: Use an appropriate amount on
a toothbrush.
Dimension: 55x179mm



Brazil's finest Propolis

Total care toothpaste containing Propolis from Brazil and Natural Ingredients for comprehensive oral care

* Propolis is a natural material that helps improve health such as immunity enhancement and antioxidant, and in the case of Brazilian products, it is receiving the best evaluation for its management on a stricter standard.



Safe toothpaste without harmful ingredients.

As it uses only **safe, natural-derived surfactants**, it might produce less foam, but it has nothing to do with cleaning power



Prevention of tooth decay

Fluorine ions (792 ppm fluoride content) form
a strong fluorine protection film to prevent tooth decay



Propolis extract



Chamomile extract



Bottlenose extract



Royal Jelly extract



Green tea extract



Aquapick Toothfoam

AQ-30 | AQ-31



Ministry of Food and
Drug Safety

Officially registered product

- Certified product **MFDS**
- Brushing & gargling at once!
- Whitening & Strengthening teeth



It can be used with
Aquapick oral irrigator.

Fill up the Aquapick water tank with
water and pump Aquapick Toothfoam
2~3 times inside. Use it after diluting
the foam with water.



Product name: Aquapick Toothfoam Model: AQ-30/AQ-31 Volume: 50ml/150ml

Smooth bubbling foam cleaning the teeth!
Innovative Toothfoam

Fine bubbles

Xylitol for prevention of dental cavities

Hydroxyapatite for white and healthy teeth

(Hydroxyapatite: a major component and
an essential ingredient of bone and teeth.
It gives bones and teeth their rigidity)

Green tea extract to
freshen up your mouth

Good for swollen and
bleeding gums
(Allatoin: prevention of
gum disease)

Contains lemon oil for shiny teeth

Abrasive-free
Surfactant-free
It is good for
toothache since it is
Abrasive-free



Fresh
mint

Whitening
effect

Removal of
bad breath

Prevent
cavity

Prevent inflammation
of the gums

Aquapick Cleaning Tablet

AQ-10



Same as actual size




Product Name: Aquapick Cleaning Tablet
Model: AQ-10
Use: For cleaning oral irrigator
Recommended frequency of use: More than 2~4 times per month
Storage: Airtight container, room temperature (0~30℃)
How to use: Put 1 tablet in 200~600ml of water and rinse 1~2 times
Amount: 12ea/package

● Cleaning even invisible spots!



✓ Fur in water tank

✓ Water streaming spots



✓ Substance within the hose

✓ Invisible spots



● Great sterilizing effect

Sterilize 99.9% of Escherichia coli, Staphylococcus and Mutans
- Certificate authority: Korea Testing & Research Institute (KTR)



● Made safe

Alcohol-free, Chlorine-free, Dye-free
Individual packaging facilitates easy use



● Sanitary cleaning

Even after cleaning once, it gets as clean as drinking water
- Certificate authority: Korea Environment & Water Works Institute (KEWWI)



● Safe Type 2 Cleaner

Type 2 cleaner is used for kitchenware cleaning
- Certificate authority: Korea Conformity Laboratories (KCL)

● With one Aquapick cleaning tablet, it takes care of the internal cleaning for oral irrigator!



▶ Oral irrigator for home use



▶ Portable oral irrigator



“Improving Oral Health Globally,”

Aquapick has achieved international **approval for safety and quality**



EU certification



Korean certification
for safe electric
appliances



Japanese
certification



European certification for
meeting the approved
electrowave standard



Recommended by
the Korean academy
of oral health



Korean intellectual
Property office